

## Product datasheet for **GA115706**

### Human TRIM67-AS1 activation kit by CRISPRa

#### Product data:

|                      |                                                                                                                                                                                                                                                                                                                           |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b> | CRISPRa Kits                                                                                                                                                                                                                                                                                                              |
| <b>Description:</b>  | LOC149373 CRISPRa kit - CRISPR gene activation of human uncharacterized LOC149373                                                                                                                                                                                                                                         |
| <b>Symbol:</b>       | TRIM67-AS1                                                                                                                                                                                                                                                                                                                |
| <b>Components:</b>   | GA115706G1, TRIM67-AS1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA115706G2, TRIM67-AS1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA115706G3, TRIM67-AS1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| <b>Format:</b>       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                                   |

#### Product images:

