

Product datasheet for **GA115708**

Human BNIP1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	BNIP1 CRISPRa kit - CRISPR gene activation of human BCL2 interacting protein like
Symbol:	BNIP1
Components:	GA115708G1, BNIP1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115708G2, BNIP1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115708G3, BNIP1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q7Z465</u>
Synonyms:	BNIP-S; BNIP-Salpha; BNIP-Sbeta; BNIP1-1; BNIP1-2; BNIP1; BNIP2; BNIPS; PP753
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)