

Product datasheet for **GA115716**

Human SRARP activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SRARP CRISPRa kit - CRISPR gene activation of human steroid receptor associated and regulated protein
Symbol:	SRARP
Components:	GA115716G1, SRARP gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115716G2, SRARP gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115716G3, SRARP gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8NEQ6</u>
Synonyms:	ERRF
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)