

Product datasheet for **GA115732**

Human Snf1k (SIK1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SIK1 CRISPRa kit - CRISPR gene activation of human salt inducible kinase 1
Symbol:	SIK1
UniProt ID:	P57059
Synonyms:	MSK; SIK; SNF1LK
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)