

Product datasheet for **GA115735**

Human **XKR3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	XKR3 CRISPRa kit - CRISPR gene activation of human XK related 3
Symbol:	XKR3
Components:	GA115735G1, XKR3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115735G2, XKR3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115735G3, XKR3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5GH77</u>
Synonyms:	XRG3; XTES
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

