

Product datasheet for **GA115743**

Human DUSP18 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DUSP18 CRISPRa kit - CRISPR gene activation of human dual specificity phosphatase 18
Symbol:	DUSP18
UniProt ID:	<u>Q8NEJ0</u>
Synonyms:	DSP18; DUSP20; LMWDSP20
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)