


Product datasheet for **GA115749**

Human PHETA2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PHETA2 CRISPRa kit - CRISPR gene activation of human PH domain containing endocytic trafficking adaptor 2
Symbol:	PHETA2
Components:	GA115749G1, PHETA2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115749G2, PHETA2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115749G3, PHETA2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q6ICB4
Synonyms:	IPIP27B; Ses2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)