

Product datasheet for **GA115773**

Human PUS10 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PUS10 CRISPRa kit - CRISPR gene activation of human pseudouridine synthase 10
Symbol:	PUS10
Components:	GA115773G1, PUS10 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115773G2, PUS10 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115773G3, PUS10 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3MIT2</u>
Synonyms:	CCDC139; DOBI
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)