

## Product datasheet for **GA115779**

### Human ERFE activation kit by CRISPRa

#### Product data:

|                      |                                                                                                                                                                                                                                                                                                         |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b> | CRISPRa Kits                                                                                                                                                                                                                                                                                            |
| <b>Description:</b>  | ERFE CRISPRa kit – CRISPR gene activation of human erythroferrone                                                                                                                                                                                                                                       |
| <b>Symbol:</b>       | ERFE                                                                                                                                                                                                                                                                                                    |
| <b>Components:</b>   | GA115779G1, ERFE gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA115779G2, ERFE gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA115779G3, ERFE gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| <b>UniProt ID:</b>   | <u><a href="#">Q4G0M1</a></u>                                                                                                                                                                                                                                                                           |
| <b>Synonyms:</b>     | CIQTNF15; CTRP15; FAM132B                                                                                                                                                                                                                                                                               |
| <b>Format:</b>       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                 |

#### Product images:

