


Product datasheet for **GA115787**

Human SLC38A11 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SLC38A11 CRISPRa kit - CRISPR gene activation of human solute carrier family 38 member 11
Symbol:	SLC38A11
Components:	GA115787G1, SLC38A11 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115787G2, SLC38A11 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115787G3, SLC38A11 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q08AI6</u>
Synonyms:	AVT2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)