

Product datasheet for **GA115788**

Human SLC23A3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SLC23A3 CRISPRa kit – CRISPR gene activation of human solute carrier family 23 member 3
Symbol:	SLC23A3
Components:	GA115788G1, SLC23A3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115788G2, SLC23A3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115788G3, SLC23A3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6PIS1</u>
Synonyms:	E2BP3; SVCT3; Yspll
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

