

Product datasheet for **GA115801**

Human TFAA4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TFAA4 CRISPRa kit - CRISPR gene activation of human TFAA chemokine like family member 4
Symbol:	TFAA4
UniProt ID:	Q96LR4
Synonyms:	TFAA-4; TFAA4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)