

Product datasheet for **GA115803**

Human PP2D1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PP2D1 CRISPRa kit – CRISPR gene activation of human protein phosphatase 2C like domain containing 1
Symbol:	PP2D1
Components:	GA115803G1, PP2D1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115803G2, PP2D1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115803G3, PP2D1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A8MPX8</u>
Synonyms:	C3orf48
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

