

## Product datasheet for **GA115852**

### Human DAB2IP activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | DAB2IP CRISPRa kit - CRISPR gene activation of human DAB2 interacting protein |
| Symbol:       | DAB2IP  |
| UniProt ID:   | <u><a href="#">Q5VWQ8</a></u>   |
| Synonyms:     | AF9Q34; AIP-1; AIP1; DIP1/2   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)       |

#### Product images:



[View online >](#)