

Product datasheet for **GA115854**

Human SLC36A2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SLC36A2 CRISPRa kit - CRISPR gene activation of human solute carrier family 36 member 2
Symbol:	SLC36A2
Components:	GA115854G1, SLC36A2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115854G2, SLC36A2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115854G3, SLC36A2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q495M3
Synonyms:	PAT2; TRAMD1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)