

Product datasheet for **GA115859**

Human TMEM167A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM167A CRISPRa kit - CRISPR gene activation of human transmembrane protein 167A
Symbol:	TMEM167A
Components:	GA115859G1, TMEM167A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115859G2, TMEM167A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115859G3, TMEM167A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8TBQ9</u>
Synonyms:	TMEM167
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)