

## Product datasheet for **GA115863**

### Human PLEKHG4B activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                     |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                        |
| Description:  | PLEKHG4B CRISPRa kit - CRISPR gene activation of human pleckstrin homology and RhoGEF domain containing G4B                                                                                                                                                                                                         |
| Symbol:       | PLEKHG4B                                                                                                                                                                                                                                                                                                            |
| Components:   | GA115863G1, PLEKHG4B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA115863G2, PLEKHG4B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA115863G3, PLEKHG4B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q96PX9</a></u>                                                                                                                                                                                                                                                                                       |
| Synonyms:     | FLJ43173; KIAA1909                                                                                                                                                                                                                                                                                                  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                             |

#### Product images:



[View online »](#)