


Product datasheet for **GA115871**

Human TTC23L activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TTC23L CRISPRa kit - CRISPR gene activation of human tetratricopeptide repeat domain 23 like
Symbol:	TTC23L
Components:	GA115871G1, TTC23L gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115871G2, TTC23L gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115871G3, TTC23L gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6PF05</u>
Synonyms:	MC25-1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)