

Product datasheet for **GA115876**

Human RNF158 (SH3RF2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SH3RF2 CRISPRa kit - CRISPR gene activation of human SH3 domain containing ring finger 2
Symbol:	SH3RF2
UniProt ID:	<u>Q8TEC5</u>
Synonyms:	HEPP1; POSHER; PPP1R39; RNF158
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)