

Product datasheet for **GA115878**

Human RNF145 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	RNF145 CRISPRa kit - CRISPR gene activation of human ring finger protein 145
Symbol:	RNF145
Components:	GA115878G1, RNF145 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115878G2, RNF145 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115878G3, RNF145 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q96MT1
Synonyms:	DKFZp686M11215; FLJ31951; FLJ44310
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)