

Product datasheet for **GA115881**

Human CNKSR3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CNKSR3 CRISPRa kit - CRISPR gene activation of human CNKSR family member 3
Symbol:	CNKSR3
UniProt ID:	Q6P9H4
Synonyms:	MAGI1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)