

Product datasheet for **GA115883**

Human SAMD3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SAMD3 CRISPRa kit - CRISPR gene activation of human sterile alpha motif domain containing 3
Symbol:	SAMD3
Components:	GA115883G1, SAMD3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115883G2, SAMD3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115883G3, SAMD3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8N6K7</u>
Synonyms:	FLJ34563; MGC35163
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)