

Product datasheet for **GA115890**

Human Ecat1 (KHDC3L) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KHDC3L CRISPRa kit - CRISPR gene activation of human KH domain containing 3 like, subcortical maternal complex member
Symbol:	KHDC3L
Components:	GA115890G1, Ecat1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115890G2, Ecat1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115890G3, Ecat1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q587J8</u>
Synonyms:	C6orf221; ECAT1; HYDM2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

