

Product datasheet for **GA115905**

Human TMEM213 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM213 CRISPRa kit - CRISPR gene activation of human transmembrane protein 213
Symbol:	TMEM213
Components:	GA115905G1, TMEM213 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115905G2, TMEM213 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115905G3, TMEM213 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A2RRL7</u>
Synonyms:	transmembrane protein 213
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)