

Product datasheet for **GA115917**

Human PRAGMIN (PRAG1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PRAG1 CRISPRa kit - CRISPR gene activation of human PEAK1 related, kinase-activating pseudokinase 1
Symbol:	PRAG1
UniProt ID:	Q86YV5
Synonyms:	PRAGMIN
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)