


Product datasheet for **GA115918**

Human PEBP4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PEBP4 CRISPRa kit - CRISPR gene activation of human phosphatidylethanolamine binding protein 4
Symbol:	PEBP4
Components:	GA115918G1, PEBP4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115918G2, PEBP4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115918G3, PEBP4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96S96</u>
Synonyms:	CORK-1; CORK1; GWTM1933; HEL-S-300; hPEBP4; PEBP-4; PRO4408
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)