


## Product datasheet for **GA115923**

### Human ESCO2 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                            |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                               |
| Description:  | ESCO2 CRISPRa kit - CRISPR gene activation of human establishment of sister chromatid cohesion N-acetyltransferase 2                                                                                                                                                                                       |
| Symbol:       | ESCO2                                                                                                                                                                                                                                                                                                      |
| Components:   | GA115923G1, ESCO2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA115923G2, ESCO2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA115923G3, ESCO2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q56NI9</a></u>                                                                                                                                                                                                                                                                              |
| Synonyms:     | 2410004I17Rik; EFO2; RBS                                                                                                                                                                                                                                                                                   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                    |

#### Product images:



[View online »](#)