

Product datasheet for **GA115931**

Human TMEM74 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM74 CRISPRa kit - CRISPR gene activation of human transmembrane protein 74
Symbol:	TMEM74
Components:	GA115931G1, TMEM74 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115931G2, TMEM74 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115931G3, TMEM74 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96NL1</u>
Synonyms:	NET36
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)