

Product datasheet for **GA115933**

Human C8orf48 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	C8orf48 CRISPRa kit - CRISPR gene activation of human chromosome 8 open reading frame 48
Symbol:	C8orf48
Components:	GA115933G1, C8orf48 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115933G2, C8orf48 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115933G3, C8orf48 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96LL4</u>
Synonyms:	FLJ25402
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)