

## Product datasheet for **GA115939**

### Human CAMSAP1 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	CAMSAP1 CRISPRa kit - CRISPR gene activation of human calmodulin regulated spectrin associated protein 1
Symbol:	CAMSAP1
UniProt ID:	<a href="#">Q5T5Y3</a>
Synonyms:	bA100C15.1; PRO2405
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)