

Product datasheet for **GA115949**

Human RASEF activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	RASEF CRISPRa kit - CRISPR gene activation of human RAS and EF-hand domain containing
Symbol:	RASEF
Components:	GA115949G1, RASEF gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115949G2, RASEF gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115949G3, RASEF gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8IZ41</u>
Synonyms:	RAB45
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)