


Product datasheet for **GA115977**

Human DGAT2L3 (AWAT1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	AWAT1 CRISPRa kit - CRISPR gene activation of human acyl-CoA wax alcohol acyltransferase 1
Symbol:	AWAT1
Components:	GA115977G1, DGAT2L3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115977G2, DGAT2L3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115977G3, DGAT2L3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q58HT5
Synonyms:	DGA2; DGAT2L3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)