

Product datasheet for **GA115983**

Human CXorf38 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CXorf38 CRISPRa kit - CRISPR gene activation of human chromosome X open reading frame 38
Symbol:	CXorf38
Components:	GA115983G1, CXorf38 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115983G2, CXorf38 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115983G3, CXorf38 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8TB03</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)