

Product datasheet for **GA115990**

Human **SLC35G1** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SLC35G1 CRISPRa kit - CRISPR gene activation of human solute carrier family 35 member G1
Symbol:	SLC35G1
Components:	GA115990G1, SLC35G1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115990G2, SLC35G1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115990G3, SLC35G1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q2M3R5
Synonyms:	C10orf60; TMEM20
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)