

Product datasheet for **GA116007**

Human **SLC5A8** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SLC5A8 CRISPRa kit - CRISPR gene activation of human solute carrier family 5 member 8
Symbol:	SLC5A8
Components:	GA116007G1, SLC5A8 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116007G2, SLC5A8 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116007G3, SLC5A8 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8N695</u>
Synonyms:	AIT; SMCT; SMCT1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)