

Product datasheet for **GA116008**

Human PPTC7 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | PPTC7 CRISPRa kit - CRISPR gene activation of human PPTC7 protein phosphatase homolog |
| Symbol: | PPTC7 |
| UniProt ID: | Q8NI37 |
| Synonyms: | TA-PP2C; TAPP2C |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)