

Product datasheet for **GA116016**

Human TMEM229B activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM229B CRISPRa kit - CRISPR gene activation of human transmembrane protein 229B
Symbol:	TMEM229B
Components:	GA116016G1, TMEM229B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116016G2, TMEM229B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116016G3, TMEM229B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8NBD8</u>
Synonyms:	C14orf83
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)