

Product datasheet for **GA116017**

Human Nesprin3 (SYNE3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SYNE3 CRISPRa kit - CRISPR gene activation of human spectrin repeat containing nuclear envelope family member 3
Symbol:	SYNE3
UniProt ID:	Q6ZMZ3
Synonyms:	C14orf49; Nesp3; NET53
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)