

## Product datasheet for **GA116047**

### Human RETREG3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	RETREG3 CRISPRa kit - CRISPR gene activation of human reticulophagy regulator family member 3
Symbol:	RETREG3
Components:	GA116047G1, RETREG3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116047G2, RETREG3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116047G3, RETREG3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q86VR2</u>
Synonyms:	FLJ33806
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)