

Product datasheet for **GA116051**

Human TRPV3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TRPV3 CRISPRa kit - CRISPR gene activation of human transient receptor potential cation channel subfamily V member 3
Symbol:	TRPV3
UniProt ID:	<u>Q8NET8</u>
Synonyms:	FNEPPK2; OLMS; VRL3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)