


Product datasheet for **GA116063**

Human ZNF497 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF497 CRISPRa kit - CRISPR gene activation of human zinc finger protein 497
Symbol:	ZNF497
Components:	GA116063G1, ZNF497 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116063G2, ZNF497 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116063G3, ZNF497 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6ZNH5</u>
Synonyms:	FLJ42722; FLJ44773
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)