

Product datasheet for **GA116070**

Human ZNF791 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF791 CRISPRa kit – CRISPR gene activation of human zinc finger protein 791
Symbol:	ZNF791
Components:	GA116070G1, ZNF791 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116070G2, ZNF791 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116070G3, ZNF791 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q3KP31
Synonyms:	FLJ90396; MGC126179; MGC126180
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

