

Product datasheet for **GA116076**

Human ZNF383 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF383 CRISPRa kit – CRISPR gene activation of human zinc finger protein 383
Symbol:	ZNF383
Components:	GA116076G1, ZNF383 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116076G2, ZNF383 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116076G3, ZNF383 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8NA42</u>
Synonyms:	HSD17
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

