

Product datasheet for **GA116082**

Human C19orf46 (SYNE4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SYNE4 CRISPRa kit - CRISPR gene activation of human spectrin repeat containing nuclear envelope family member 4
Symbol:	SYNE4
UniProt ID:	Q8N205
Synonyms:	C19orf46; DFNB76; Nesp4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)