

Product datasheet for **GA116083**

Human ZNF676 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF676 CRISPRa kit – CRISPR gene activation of human zinc finger protein 676
Symbol:	ZNF676
Components:	GA116083G1, ZNF676 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116083G2, ZNF676 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116083G3, ZNF676 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8N7Q3</u>
Synonyms:	zinc finger protein 676
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

