

Product datasheet for **GA116088**

Human PAP2D (PLPPR5) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLPPR5 CRISPRa kit - CRISPR gene activation of human phospholipid phosphatase related 5
Symbol:	PLPPR5
Components:	GA116088G1, PAP2D gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116088G2, PAP2D gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116088G3, PAP2D gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q32ZL2</u>
Synonyms:	PAP2; PAP2D; PRG5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)