

Product datasheet for **GA116095**

Human **CITED4** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	CITED4 CRISPRa kit - CRISPR gene activation of human Cbp/p300 interacting transactivator with Glu/Asp rich carboxy-terminal domain 4
Symbol:	CITED4
UniProt ID:	Q96RK1
Synonyms:	Cbp/p300-interacting transactivator, with Glu/Asp-rich carboxy-terminal domain, 4; OTTHUMP0000009206; transcriptional co-activator 4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)