

Product datasheet for **GA116104**

Human MPRA (PAQR7) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PAQR7 CRISPRa kit - CRISPR gene activation of human progesterin and adipoQ receptor family member 7
Symbol:	PAQR7
Components:	GA116104G1, MPRA gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116104G2, MPRA gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116104G3, MPRA gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q86WK9</u>
Synonyms:	MPRA; mSR; PGLP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)