

Product datasheet for **GA116105**

Human **TTC24** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	TTC24 CRISPRa kit - CRISPR gene activation of human tetratricopeptide repeat domain 24
Symbol:	TTC24
Components:	GA116105G1, TTC24 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116105G2, TTC24 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116105G3, TTC24 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A2A3L6</u>
Synonyms:	FLJ20249
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)