

Product datasheet for **GA116110**

Human LRRN4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LRRN4 CRISPRa kit - CRISPR gene activation of human leucine rich repeat neuronal 4
Symbol:	LRRN4
UniProt ID:	<u>Q8WUT4</u>
Synonyms:	C20orf75; dj1056H1.1; NLRR-4; NLRR4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)